

GFAAS Sample Preparation
CHEM 322
Fall 2009

Protocol:

All glass and plastic-ware used for solution preparation was cleaned by washing with dilute nitric acid. Solutions were prepared by digesting an accurately weighed sample of either NIST Bovine muscle (RM # 8414) or NIST Bovine Liver (RM # 1577b). Briefly, the digestion procedure involved addition of a 7 mL of trace metal grade concentrated nitric acid and 2 mL of 18 MΩ water to the pre-weighed samples in a Teflon microwave vessel. The vessels were inserted into the CEM MDS2000 and run through a heating cycle to promote decomposition. Two reagent blanks were also prepared following this protocol. After cooling, the samples were quantitatively transferred to clean, pre-weighed plastic bottles.

Sample Information:

ID	Mass of Sample	Mass of Empty Bottle	Mass of Bottle + Solution
Sample 1	0.6210 g	12.5030 g	44.2489 g
Sample 2	0.6688 g	12.4563 g	42.4812 g
Sample 3	1.5619 g	12.3857 g	45.1804 g
Sample 4	1.5203 g	12.4976 g	44.3273 g
Blank 1	---	12.5140 g	43.5057 g
Blank 2	---	12.4918 g	41.7392 g

Characteristics of Bovine Standards:

The Certificate of Analysis for the Bovine Muscle and Bovine Liver Standard Reference Materials can be found online at the NIST (National Institute of Standards and Technology) (www.nist.gov). Begin by clicking the "Standard Reference Materials" link. On the page that opens, click on "Search the Online Catalog". Search for the corresponding reference material using its RM# and look for either a "Certificate of Analysis" or "Report of Investigation".